



THE PATENTED

GROUT SAMPLE BOX

Proven by tests to yield comparable compressive strength values to samples cast by traditional method.


CAT#: GSB

25 per box

THE GROUT SAMPLE BOX provides specimens that are ready for testing. It yields flat surfaces with perpendicular sides to the tops and bottoms, whereas the pinwheel method usually yields specimens that need to be trimmed and cut before testing. Meets ASTM C1019, Section 6.



New interlocking tabs offer foolproof assembly and ensures a truly square box.



New interlocking tabs ensure tight fitting corners for a precise specimen.



Staggered folds are integrated to give our sample box a truly flat bottom.



New single perforated line makes it easier to fold and assemble.



