




# BUCKET INSERT


**THE BUCKET INSERT** permits a standard 3-1/2 or 5 gallon bucket to carry up to three 4x8" specimens in an upright position and separated to prevent specimen damage.

- An easy and inexpensive way to convert a standard bucket into a convenient test specimen carryall!
- Sold in cartons, 4 pieces or 12 pieces per carton
- Reusable
- Reusable bucket not included but are available at most home improvement and hardware stores



**CAT#: TC4-BI**

**CAT#: TC4-BI-4**



Patent Applied



SCAN ME



[www.deslinc.com](http://www.deslinc.com)



800-743-4106



/DeslauriersInc