



# LIFTING HANDLE

- Allows you to carry 4" and 6" specimens with one hand.
- Use to remove 4" and 6" specimens from curing tank.
- Safe and reduce back stress.
- Lightweight, easy to transport and safe, no sharp edges or corners.



**CAT#: TC6LH** --For 6" specimens  
**CAT#: TC4LH** --For 4" specimens



# 6CH ECONO CARRY HANDLE

The handy 6CH Econo Cylinder Carrier will fit any 6" diameter plastic cylinder mold with a wide lip and can be used with lids in its place. Flexible plastic ring slides onto mold from bottom prior to use and is held in place by top lip. When ready to transport, carrier handles are rotated up to snap together in the middle. Recommended for single use. Packaged in cartons of 50.

**CAT#: ECH**



cylinder NOT included



 [www.deslinc.com](http://www.deslinc.com)  
Web

 800-743-4106  
Phone

 /DeslauriersInc  
Facebook